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# HOW TO USE WEIGHTS CORRECTLY!

## Store weights in original box

All weights must come in correct aluminum box or plastic box. Box will keep your weights from gathering dust. And different sizes will help prevent mix-ups.

## Clean before using

Use brush and wipes to remove rust, which can confirm the accuracy for weights.

## Put weights beside the balance

Putting weights beside the balance will help to ensure that balance and weights are similarly acclimatized.

## Use tweezers for small weights.

Use tweezers with coated tips to lift 1mg-200g weights.

## Inspect for damage before using

Before using the weights, check the surface for impurity, corrosion or water. Pay more attention to dust, fingerprints and surface marks.

## Use handles for big weights.

Use specialized handles for big weights not only protects the weight's surface finish, but also provides a more ergonomic lifting experience to protect the health of the operator.

## Do not touch weights with hands

This is especially critical for weights with smaller nominal weights. Skin acids disrupt weight's surface.

## Do not forget periodic recalibration

Adhering to recalibration schedules stipulated by ISO9001 will help keep process accurate. Weights need to be re-calibrated each year.